

Product Datasheet

Mouse bACE1 ELISA Kit (Ready to Use) (orb553036)

Description	Mouse Beta-Site APP Cleaving Enzyme 1 (bACE1) Sandwich ELISA KitReady to Use: Pre-diluted detection ...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cerebrospinal fluid or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cerebrospinal fluid and other biological fluids.
Uniprot ID	P56818
Sensitivity	5.9pg/mL
Entrez	23821
Expiration Date	Please enquire.