

Product Datasheet

Mouse CACNa1S ELISA Kit (Ready to Use) (orb553073)

Description	Mouse Calcium Channel, Voltage Dependent, L-Type, Alpha 1S Subunit (CACNa1S) Sandwich ELISA KitReady...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q02789
Sensitivity	0.101ng/mL
Entrez	12292
Expiration Date	Please enquire.