

Product Datasheet

Mouse CT1 ELISA Kit (Ready to Use) (orb553094)

Description	Mouse Cardiotrophin 1 (CT1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conve...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q60753
Sensitivity	6.2pg/mL
Entrez	13019
Expiration Date	Please enquire.