

Product Datasheet

Mouse CLCF1 ELISA Kit (Ready to Use) (orb553095)

Description	Mouse Cardiotrophin Like Cytokine Factor 1 (CLCF1) Sandwich ELISA KitReady to Use: Pre-diluted detec...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9QZM3
Sensitivity	6.4pg/mL
Entrez	56708
Expiration Date	Please enquire.