

Product Datasheet

Mouse Caspase 8 ELISA Kit (Ready to Use) (orb553105)

Description	Mouse Caspase 8 (Caspase 8) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for conve...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	089110
Sensitivity	0.069ng/mL
Entrez	12370
Expiration Date	Please enquire.