

Product Datasheet

Mouse CTSB ELISA Kit (Ready to Use) (orb553112)

Description	Mouse Cathepsin B (CTSB) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenie...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P10605
Sensitivity	0.057ng/mL
Entrez	13030
Expiration Date	Please enquire.