

Product Datasheet

Mouse CCKAR ELISA Kit (Ready to Use) (orb553147)

Description	Mouse Cholecystokinin A Receptor (CCKAR) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reage...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	O08786
Sensitivity	0.064ng/mL
Entrez	12425
Expiration Date	Please enquire.