

Product Datasheet

Mouse CPP ELISA Kit (Ready to Use) (orb553215)

Description	Mouse Copeptin (CPP) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient a...
Reactivity	Mouse
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P35455
Sensitivity	8.83pg/mL
Entrez	11998
Expiration Date	Please enquire.