

Product Datasheet

Mouse CRN ELISA Kit (Ready to Use) (orb553216)

Description	Mouse Corin (CRN) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient and ...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9Z319
Sensitivity	6.7pg/mL
Entrez	53419
Expiration Date	Please enquire.