

## Product Datasheet

### Mouse NTXI ELISA Kit (Ready to Use) (orb553224)

<b>Description</b>	Mouse Cross Linked N-Telopeptide Of Type I Collagen (NTXI) Sandwich ELISA KitReady to Use: Pre-dilut...
<b>Reactivity</b>	Mouse
<b>Range</b>	3.12-200ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilityserum, plasma, urine or other biological fluids.
<b>Sample Types</b>	serum, plasma, urine and other biological fluids.
<b>Sensitivity</b>	0.96ng/mL
<b>Expiration Date</b>	Please enquire.