

Product Datasheet

Mouse CYP27B1 ELISA Kit (Ready to Use) (orb553242)

Description	Mouse Cytochrome P450 27B1 (CYP27B1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents ...
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O35084
Sensitivity	0.214ng/mL
Entrez	13115
Expiration Date	Please enquire.