

Product Datasheet

Mouse DKK3 ELISA Kit (Ready to Use) (orb553268)

Description	Mouse Dickkopf Related Protein 3 (DKK3) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagen...
Reactivity	Mouse
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	14.3pg/mL
Expiration Date	Please enquire.