

## Product Datasheet

### Mouse EGFL7 ELISA Kit (Ready to Use) (orb553298)

<b>Description</b>	Mouse EGF Like Domain Protein, Multiple 7 (EGFL7) Sandwich ELISA KitReady to Use: Pre-diluted detect...
<b>Reactivity</b>	Mouse
<b>Range</b>	31.2-2000pg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9QXT5</a>
<b>Sensitivity</b>	15.5pg/mL
<b>Entrez</b>	<a href="#">353156</a>
<b>Expiration Date</b>	Please enquire.