

Product Datasheet

Mouse Neutrophil Elastase ELISA Kit (Ready to Use) (orb553301)

Description	Mouse Elastase 2, Neutrophil (Neutrophil Elastase) Sandwich ELISA KitReady to Use: Pre-diluted detec...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q3UP87
Sensitivity	0.063ng/mL
Entrez	50701
Expiration Date	Please enquire.