

Product Datasheet

Mouse FABP6 ELISA Kit (Ready to Use) (orb553342)

| | |
|--------------------------|--|
| Description | Mouse Fatty Acid Binding Protein 6, Ileal (FABP6) Sandwich ELISA KitReady to Use: Pre-diluted detect... |
| Reactivity | Mouse |
| Range | 78.1-5000pg/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids. |
| Uniprot ID | P51162 |
| Sensitivity | 23pg/mL |
| Entrez | 16204 |
| Expiration Date | Please enquire. |