

Product Datasheet

Mouse FGF15 ELISA Kit (Ready to Use) (orb553361)

| | |
|--------------------------|--|
| Description | Mouse Fibroblast Growth Factor 15 (FGF15) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagent... |
| Reactivity | Mouse |
| Range | 31.2-2000pg/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| Uniprot ID | O35622 |
| Sensitivity | 12.4pg/mL |
| Entrez | 14170 |
| Expiration Date | Please enquire. |