


Product Datasheet

Mouse FSH ELISA Kit (Ready to Use) (orb553382)

Description	Mouse Follicle Stimulating Hormone (FSH) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagent...
Reactivity	Mouse
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q60687
Sensitivity	1.49ng/mL
Entrez	14308
Expiration Date	Please enquire.



Standard Curve of Mouse FSH ELISA Kit.