

Product Datasheet

Mouse GPER ELISA Kit (Ready to Use) (orb553394)

Description	Mouse G Protein Coupled Estrogen Receptor 1 (GPER) Sandwich ELISA KitReady to Use: Pre-diluted detec...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q8BMP4
Sensitivity	0.054ng/mL
Entrez	76854
Expiration Date	Please enquire.