

Product Datasheet

Mouse GLS ELISA Kit (Ready to Use) (orb553437)

Description	Mouse Glutaminase (GLS) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenien...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	D3Z7P3
Sensitivity	0.117ng/mL
Entrez	14660
Expiration Date	Please enquire.