

Product Datasheet

Mouse GLS2 ELISA Kit (Ready to Use) (orb553438)

Description	Mouse Glutaminase 2 (GLS2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q571F8
Sensitivity	0.25ng/mL
Entrez	216456
Expiration Date	Please enquire.