

Product Datasheet

Mouse GHRH ELISA Kit (Ready to Use) (orb553480)

Description	Mouse Growth Hormone Releasing Hormone (GHRH) Sandwich ELISA Kit Ready to Use: Pre-diluted detection ...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P16043
Sensitivity	4.34pg/mL
Entrez	14601
Expiration Date	Please enquire.