

Product Datasheet

Mouse HSPA8 ELISA Kit (Ready to Use) (orb553490)

Description	Mouse Heat Shock 70kDa Protein 8 (HSPA8) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reage...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P63017
Sensitivity	0.116ng/mL
Entrez	15481
Expiration Date	Please enquire.