

Product Datasheet

Mouse HMG1 ELISA Kit (Ready to Use) (orb553509)

| | |
|--------------------------|---|
| Description | Mouse High Mobility Group Protein 1 (HMG1) Sandwich ELISA KitReady to Use: Pre-diluted detection rea... |
| Reactivity | Mouse |
| Range | 62.5-4000pg/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Uniprot ID | P63158 |
| Sensitivity | 23.73pg/mL |
| Entrez | 15289 |
| Expiration Date | Please enquire. |