

## Product Datasheet

### Mouse HADHb ELISA Kit (Ready to Use) (orb553530)

|                          |                                                                                                          |
|--------------------------|----------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Hydroxyacyl Coenzyme A Dehydrogenase Beta (HADHb) Sandwich ELISA Kit Ready to Use: Pre-diluted ... |
| <b>Reactivity</b>        | Mouse                                                                                                    |
| <b>Range</b>             | 0.312-20ng/mL                                                                                            |
| <b>Form/Appearance</b>   | Ready to Use                                                                                             |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt                                                      |
| <b>Note</b>              | For research use only                                                                                    |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stability tissue homogenates or other biological fluids.     |
| <b>Sample Types</b>      | tissue homogenates and other biological fluids.                                                          |
| <b>Sensitivity</b>       | 0.115ng/mL                                                                                               |
| <b>Expiration Date</b>   | Please enquire.                                                                                          |