

Product Datasheet

Mouse IDE ELISA Kit (Ready to Use) (orb553552)

Description	Mouse Insulin Degrading Enzyme (IDE) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Mouse
Range	3.12-200pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9JHR7
Sensitivity	1.33pg/mL
Entrez	15925
Expiration Date	Please enquire.