

Product Datasheet

Mouse INSL6 ELISA Kit (Ready to Use) (orb553560)

Description	Mouse Insulin Like Protein 6 (INSL6) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9QY05
Sensitivity	30pg/mL
Entrez	27356
Expiration Date	Please enquire.