

Product Datasheet

Mouse ITGa1 ELISA Kit (Ready to Use) (orb553564)

Description	Mouse Integrin Alpha 1 (ITGa1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for co...
Reactivity	Mouse
Range	6.25-400pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q3V3R4
Sensitivity	2.65pg/mL
Entrez	109700
Expiration Date	Please enquire.