

Product Datasheet

Mouse ITGa3 ELISA Kit (Ready to Use) (orb553565)

Description	Mouse Integrin Alpha 3 (ITGa3) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for co...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q62470
Sensitivity	0.118ng/mL
Entrez	16400
Expiration Date	Please enquire.