

Product Datasheet

Mouse ITLN1 ELISA Kit (Ready to Use) (orb553569)

Description	Mouse Intelectin 1 (ITLN1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Mouse
Range	2.5-160pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	D9J349
Sensitivity	0.9pg/mL
Entrez	16429
Expiration Date	Please enquire.