

## Product Datasheet

# Mouse IL13Ra2 ELISA Kit (Ready to Use) (orb553592)

<b>Description</b>	Mouse Interleukin 13 Receptor Alpha 2 (IL13Ra2) Sandwich ELISA KitReady to Use: Pre-diluted detectio...
<b>Reactivity</b>	Mouse
<b>Range</b>	0.312-20ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<b>O88786</b>
<b>Sensitivity</b>	0.126ng/mL
<b>Entrez</b>	<b>16165</b>
<b>Expiration Date</b>	Please enquire.