

Product Datasheet

Mouse IL5Ra ELISA Kit (Ready to Use) (orb553617)

Description	Mouse Interleukin 5 Receptor Alpha (IL5Ra) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Mouse
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P21183
Sensitivity	28.5pg/mL
Entrez	16192
Expiration Date	Please enquire.