

Product Datasheet

Mouse KLK1 ELISA Kit (Ready to Use) (orb553630)

Description	Mouse Kallikrein 1 (KLK1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Mouse
Range	18.7-1200pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P15947
Sensitivity	7.52pg/mL
Entrez	16612
Expiration Date	Please enquire.