

Product Datasheet

Mouse KLK4 ELISA Kit (Ready to Use) (orb553633)

Description	Mouse Kallikrein 4 (KLK4) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9Z0M1
Sensitivity	0.065ng/mL
Entrez	56640
Expiration Date	Please enquire.