

Product Datasheet

Mouse KISS1 ELISA Kit (Ready to Use) (orb553646)

Description	Mouse Kisspeptin 1 (KISS1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q6Y4S4
Sensitivity	4.46pg/mL
Entrez	280287
Expiration Date	Please enquire.