

Product Datasheet

Mouse LDHB ELISA Kit (Ready to Use) (orb553653)

Description	Mouse Lactate Dehydrogenase B (LDHB) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents ...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P16125
Sensitivity	0.065ng/mL
Entrez	16832
Expiration Date	Please enquire.