

Product Datasheet

Mouse LTF ELISA Kit (Ready to Use) (orb553654)

Description	Mouse Lactoferrin (LTF) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenien...
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, breast milk or other biological fluids.
Sample Types	serum, plasma, breast milk and other biological fluids.
Uniprot ID	P08071
Sensitivity	0.241ng/mL
Entrez	17002
Expiration Date	Please enquire.