

Product Datasheet

Mouse MIP3a ELISA Kit (Ready to Use) (orb553709)

Description	Mouse Macrophage Inflammatory Protein 3 Alpha (MIP3a) Sandwich ELISA Kit Ready to Use: Pre-diluted de...
Reactivity	Mouse
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	089093
Sensitivity	27pg/mL
Entrez	20297
Expiration Date	Please enquire.