

Product Datasheet

Mouse MIP3 beta ELISA Kit (Ready to Use) (orb553710)

Description	Mouse Macrophage Inflammatory Protein 3 Beta (MIP3 beta) Sandwich ELISA KitReady to Use: Pre-diluted...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O70460
Sensitivity	5.9pg/mL
Entrez	24047
Expiration Date	Please enquire.