

Product Datasheet

Mouse MMP15 ELISA Kit (Ready to Use) (orb553728)

| | |
|--------------------------|---|
| Description | Mouse Matrix Metalloproteinase 15 (MMP15) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagent... |
| Reactivity | Mouse |
| Range | 78.1-5000pg/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids. |
| Uniprot ID | 054732 |
| Sensitivity | 28.4pg/mL |
| Entrez | 17388 |
| Expiration Date | Please enquire. |