

## Product Datasheet

# Mouse NT-ProANP ELISA Kit (Ready to Use) (orb553796)

<b>Description</b>	Mouse N-Terminal Pro Atrial Natriuretic Peptide (NT-ProANP) Sandwich ELISA Kit Ready to Use: Pre-dilu...
<b>Reactivity</b>	Mouse
<b>Range</b>	78.1-5000pg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Sensitivity</b>	29pg/mL
<b>Expiration Date</b>	Please enquire.