

Product Datasheet

Mouse NMB ELISA Kit (Ready to Use) (orb553818)

Description	Mouse Neuromedin B (NMB) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenie...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	D3Z621
Sensitivity	4.75pg/mL
Entrez	68039
Expiration Date	Please enquire.