

Product Datasheet

Mouse NPFF ELISA Kit (Ready to Use) (orb553820)

Description	Mouse Neuropeptide FF (NPFF) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conv...
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityurine, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	urine, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9WVA8
Sensitivity	0.048ng/mL
Entrez	54615
Expiration Date	Please enquire.