

Product Datasheet

Mouse NPS ELISA Kit (Ready to Use) (orb553821)

Description	Mouse Neuropeptide S (NPS) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P0C0P8
Sensitivity	3.9pg/mL
Entrez	100043254
Expiration Date	Please enquire.