

Product Datasheet

Mouse PTMS ELISA Kit (Ready to Use) (orb553880)

Description	Mouse Parathymosin (PTMS) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Mouse
Range	46.8-3000ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9D0J8
Sensitivity	14.2ng/mL
Entrez	69202
Expiration Date	Please enquire.