

Product Datasheet

Mouse PTHrP ELISA Kit (Ready to Use) (orb553882)

Description	Mouse Parathyroid Hormone Related Protein (PTHrP) Sandwich ELISA Kit Ready to Use: Pre-diluted detect...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P22858
Sensitivity	5.6pg/mL
Entrez	19227
Expiration Date	Please enquire.