

Product Datasheet

Mouse PRDX1 ELISA Kit (Ready to Use) (orb553894)

Description	Mouse Peroxiredoxin 1 (PRDX1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for con...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P35700
Sensitivity	0.103ng/mL
Entrez	18477
Expiration Date	Please enquire.