

Product Datasheet

Mouse PRDX4 ELISA Kit (Ready to Use) (orb553895)

Description	Mouse Peroxiredoxin 4 (PRDX4) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O08807
Sensitivity	0.108ng/mL
Entrez	53381
Expiration Date	Please enquire.