

## Product Datasheet

### Mouse PDE1C ELISA Kit (Ready to Use) (orb553907)

<b>Description</b>	Mouse Phosphodiesterase 1C, Calmodulin Dependent (PDE1C) Sandwich ELISA KitReady to Use: Pre-diluted...
<b>Reactivity</b>	Mouse
<b>Range</b>	0.156-10ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
<b>Sample Types</b>	tissue homogenates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q64338</a>
<b>Sensitivity</b>	0.053ng/mL
<b>Entrez</b>	<a href="#">18575</a>
<b>Expiration Date</b>	Please enquire.