

Product Datasheet

Mouse PLA2G2A ELISA Kit (Ready to Use) (orb553916)

Description	Mouse Phospholipase A2, Group IIA (PLA2G2A) Sandwich ELISA KitReady to Use: Pre-diluted detection re...
Reactivity	Mouse
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P31482
Sensitivity	5.7pg/mL
Entrez	18780
Expiration Date	Please enquire.