

Product Datasheet

Mouse PGR ELISA Kit (Ready to Use) (orb553955)

Description	Mouse Progesterone Receptor (PGR) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for...
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	A6H6A5
Sensitivity	0.128ng/mL
Entrez	18667
Expiration Date	Please enquire.